



Code No.KA009

For research use only

# Heavy metal detoxification and Oxidative stress Anti Rabbit Metallothionein Monoclonal Antibody (Clone No.1A12)

Metallothinein(MT), a cysteine-rich(30%), metal-binding protein, exists in most tissues and is easily induced by many stimuli. It's multifunction roles in the body have been proposed such as a chelator to harmful heavy metals and excessive essential metals, a scavenger to various radicals and active oxygen species, and a regulator in the cell proliferation process.

This product consists of protein G purified mouse monoclonal antibody against rabbit MT-2 and can be used for immunoblotting. This antibody reacts with rabbit MT-1,2, human MT-1, mouse MT-1 and rat MT-1.

Package Size  $100 \mu g (100 \mu L/vial)$ 

Format Mouse monoclonal antibody 1.0mg/mL

Buffer 2% Block Ace as a stabilizer, containing 0.1% proclin as a bacteriostat

Storage Store below  $-20^{\circ}$ C.

Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.

Purification method The splenic lymphocytes from BALB/c mouse, immunized with rabbit

metallothionein-2, were fused to myeloma P3U1 cells. The cell line (1A12) with positive reaction was grown in ascitic fluid of BALB/c mouse, from which the

antibody was purified by Protein G affinity chromatography.

Working dilution for immunoblotting:  $0.01 \mu \text{ g/mL}$ 

Specificity This antibody reacts with rabbit MT-1,2, human MT-1, mouse MT-1 and rat MT-1

## [References]

- **1.** Yasutake A, et al. (1998): Induction by mercury compounds of brain metallothionein in rats: Hgo exposure induces long-lived brain metallothionein. *Arch Toxicol*. 72, 187-191
- **2.** Mullins JE, et al. (1998): Immunohistochemical detection of metallothionein in liver, duodenum and kidney after dietary copper-overload in rats. *Histol Histopathol*. 13, 627-633
- **3.** Kikuchi Y, et al. (1993): Induction of metallothionein in a human astrocytoma cell line by interleukin-1 and heavy metals. *FFBS Lett.* 317, 22-26
- 4. Nordberg M. (1998): Metallothionein: historical review and state of knowledge. *Talanta* 46, 243-254
- **5.** Akintola DF, et al. (1995): Development of an enzyme-linked immunosorbent assay for human metallothionein-1 in plasma and urine. *J Lab.* 126, 119-127
- 6 Nagamatsu S, et al. (2022): Sophisticated expression responses of ZNT1 and MT in response to changes in the expression of ZIPs. *Sci Rep.* 12, 7334
- 7 Ueda S, et al. (2022): Early secretory pathway-resident Zn transporter proteins contribute to cellular sphingolipid metabolism through activation of sphingomyelin phosphodiesterase 1. Am J Physiol Cell Physiol. 322, C948–C959

### Manufacturer



### Previous manufacturer



# **Kobe Research Institute**

7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-945-7075 FAX:+81-78-306-0694

URL:https://soyaku.co.jp/english/ tech-kobe@soyaku.co.jp