Code No.KG113



Research use only

Obesity and Metabolic Syndrome Related Antibody Anti Human PPARy Polyclonal Antibody


PPARs (peroxisome proliferator-activated receptors) are a family of transcription factors belonging to the nuclear hormone receptor superfamily. Widely expressed in vertebrates, PPARs play critical roles in metabolism and differentiation of a number of cell types.

The PPAR γ subtype was originally identified as a factor binding to a fatty acid specific enhancer of the aP2 gene. PPAR- γ actions are mediated by three isoforms resulting from alternative promoter selection and alternative splicing. PPAR- γ 1 is widely expressed while PPAR- γ 2 expression is restricted to adipose tissue and PPAR- γ 3 expression is restricted to adipose tissue, macrophage, and colon.

PPARy participates in adipose cell differentiation and energy storage (Ref.1).

Recently, these roles of PPAR γ have focused attention on PPAR γ as a target of the anti-diabetic thiazolidinedione class of drugs (Ref.2).

Package Size	$100\mu g(400\mu L/vial)$
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 Once thawed, store at 4 . Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was prepared from the serum of a rabbit immunized with a partial peptide representing the C-terminal domain of Human PPAR γ , and purified by peptide affinity chromatography.
Working dilution	For Western blotting: 1.0µg/ml



<u>Western blotting</u> Sample: Extracted proteins from mouse adipose tissue



Code No.KG113

Obesity and Metabolic Syndrome Related Antibody Anti Human PPARy Polyclonal Antibody

[Reference]

- Auwerx J. et al.: PPAR-gamma: a thrifty transcription factor. Nucl Recept Signal. 2003;1:e006. Nucl Recept Signal. 2003;1:e006.
- Kubota N. et al. : PPAR gamma mediates high-fat diet-induced adipocyte hypertrophy and insulin resistance. Mol Cell. 1999 Oct;4(4):597-609.

Manufacturer



7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047Telephone: +81-78-306-0295FAX:+81-78-306-0296URL:http://www.transgenic.co.jptechstaff@transgenic.co.jp