


## Anti Rat Aldehyde Oxidase Polyclonal Antibody

Package Size	25 $\mu$ g (100 $\mu$ L / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 Once thawed, store at 4 . Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with Rat Aldehyde Oxidase purified from Rat liver by antigen affinity chromatography.
Working dilution	For western blotting ; 0.5 ~ 2 $\mu$ g/mL



### Western Blotting

Sample :

Aldehyde Oxidase purified from Rat liver (0.04  $\mu$ g)

Xanthine Oxidase (0.1  $\mu$ g)

Preparation of antibodies and instruction :

Tanaka Y, Masubuchi A

Dept. of Pharmaceutics I. Tohoku Pharmaceutical University

## Anti Rat Aldehyde Oxidase Polyclonal Antibody

### 【Reference】

1. Beedham C. : Drug Metab Rev. 1985;16(1-2):119-56
2. Kisker C. et al. : Annu Rev Biochem. 1997;66:233-67
3. Wright RM. et al. : J Biol Chem. 1999 Feb 5;274(6):3878-86
4. Maia L. et al. : Arch Biochem Biophys. 2002 Apr 1;400(1):48-53
5. Beedham C. 'Molybdenum Hydroxylases.' *In Enzyme Systems that Metabolise Drugs and Other Xenobiotics*, Ioannides C. and John Wiley (ed.), New York, 2002, 147-87

#### Manufacturer



7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-306-0295      FAX:+81-78-306-0296

URL:<http://www.transgenic.co.jp>      [techstaff@transgenic.co.jp](mailto:techstaff@transgenic.co.jp)