


Anti Rat Podocalyxin Polyclonal Antibody

Package Size	25 μ g (100 μ L/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Below -20
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of Rat Podocalyxin by peptide affinity chromatography.
Working dilution	Western Blotting : 2 ~ 5 μ g/mL Immunofluorescence : 10 μ g/mL



Western blotting

Sample : 200 μ g/lane

Rat liver lysate (Negative control)

Rat kidney lysate

Rat heart lysate

Preparation of antibodies and instruction :

Dr. Takeda T. and Dr. Gejyo F. at Department of Applied Molecular Medicine, Niigata University Graduate School of Medical and Dental Sciences

Anti Rat Podocalyxin Polyclonal Antibody

[Reference]

1. Takeda T. et al. : J Clin Invest 2001 Jul;108(2):289-301
2. Kerjaschki D. et al. : Am J Pathol 1985 Mar;118(3):343-9
3. Kerjaschki D. et al. : J Cell Biol 1984;98:1591-6

Manufacturer



7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-306-0295 FAX:+81-78-306-0296

URL:<http://www.transgenic.co.jp> techstaff@transgenic.co.jp