


For research use only

Anti Human Mucosa Associated Lymphoid Tissue 1 Polyclonal Antibody

Package Size	25 μ g (100 μ L / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 Once thawed, store at 4 . Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with partial peptide of human Mucosa Associated Lymphoid Tissue 1 by peptide affinity chromatography.
Working dilution	For western blotting ; 0.5 ~ 2 μ g/mL



Western Blotting

Sample :
Human MALT1 (IPTG+) expressed in E Coli.
(IPTG-)

Preparation of antibodies and instruction :
Imai K, Chiba T
The Nippon Dental University

Anti Human Mucosa Associated Lymphoid Tissue 1 Polyclonal Antibody

【Reference】

1. Kaufman CK. et al. : J Cell Biol. 2000 May;149(5):999-1004
2. Willis TG. et al. : Cell 1999 Jan 8;96(1):35-45
3. Lucas PC. et al. : J Biol Chem. 2001 Jun 1;276(22):19012-9
4. Li Q. et al. : Genes Dev. 1999 May 15;13(10):1322-8
5. Seitz CS. et al. : Proc Natl Acad Sci USA. 1998 Mar 3;95(5):2307-12
6. Zandi E. et al. : Cell 1997 Oct 17;91(2):243-52

Manufacturer



7-1-14Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-306-0295 FAX:+81-78-306-0296

URL:<http://www.transgenic.co.jp> techstaff@transgenic.co.jp