


Anti Mouse Semaphorin 3B Polyclonal Antibody

Package Size	100 μ g (400 μ L / vial)
Format	Rabbit Polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 Once thawed, store at 4 . Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of Mouse Semaphorin 3B by peptide affinity chromatography.
Use	Immunohistochemistry, Immunoprecipitation, Western blotting
Working dilution	For Immunohistochemistry ; 15 μ g/mL For Western blotting ; 8 μ g/mL For Immunoprecipitation ; 2 μ g



Immunohistochemistry

Sample : Human placenta

Preparation of antibodies and instruction :
Dr. Suzuki, T. at Second Department of Pathology, School of
Medicine, Fukushima Medical University

Anti Mouse Semaphorin 3B Polyclonal Antibody

【Reference】

1. Castro-Rivera E. et al. : Proc Natl Acad Sci USA. 2004 Aug 3;101(31):11432-7
2. Tse C. et al. : Cancer Res. 2002 Jan 15;62(2):542-6
3. Gong S. et al. : Nature. 2003 Oct 30;425(6961):917-25

Manufacturer



7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-306-0295 FAX:+81-78-306-0296

URL:<http://www.transgenic.co.jp> techstaff@transgenic.co.jp