

Anti human Tensin2 Polyclonal Antibody

Package Size	50 μ g (200 μ L / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 Once thawed, store at 4 . Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of human Tensin2 by peptide affinity chromatography.
Working dilution	For Immunohistochemistry ; 5 μ g/mL For Western blotting ; 2 μ g



Immunohistochemistry

Sample : Mouse Kidney

Preparation of antibodies and instruction :

Dr.Agui T. and Sasaki N.

The Graduate School of Veterinary Medicine, Hokkaido University Department of Disease Control

Anti human Tensin2 Polyclonal Antibody

【Reference】

1. Lo SH. : Int J Biochem Cell Biol. 2004 Jan;36(1):31-4.
2. Lo SH. et al. : J Cell Biol. 1997 Mar 24;136(6):1349-61.
3. Calderwood DA et al. : Proc Natl Acad Sci U S A. 2003 Mar 4;100(5):2272-7.
4. Chen H. : Proc Natl Acad Sci U S A. 2002 Jan 22;99(2):733-8.
5. Lo SH. : Cancer Res. 2002 Aug 1;62(15):4217-21.

Manufacturer



7-1-14 Minatojimaminami-machi, Chuo-ku, Kobe, Japan 650-0047

Telephone: +81-78-306-0295 FAX:+81-78-306-0296

URL:<http://www.transgenic.co.jp> techstaff@transgenic.co.jp